



# V-Sorb 2800TP™

**Turbo Vacuum Pumped Surface Area & Porosity Analyzer**

**Gold APP Instruments China**

**Lead You to Particle World Better**



## Turbo Molecular Vacuum Pumped

### Surface Area & Porosity Analyzer V-SORB 2800TP Parameters

**Analysis Method:** Static volumetric nitrogen adsorption principle

**Versatility:** Adsorption and desorption isotherms, Single and multi-point BET surface area, Langmuir surface area, Average particle diameter estimation, T-plot external surface area, True density analysis, BJH mesopore total pore volume and pore size distribution, T-plot micropore, Horvath-Kawazoe (HK) micropore, Saito-Foley (SF) micropore, Dubinin-Radushkevich (DR) micropore, Dubinin-Astakhov (DA) micropore, carbon black STSA etc.

**Vacuum System:** V-Sorb unique *stainless steel micro welding vacuum manifolds system, metal-to-metal (VCR) fittings*, reduce dead space in maximum degree and also keep long duration high vacuum; the *VCR metal fittings' vacuum leak rate can reach  $1 \times 10^{-10}$  (Pa\*m<sup>3</sup>/s)*, absolutely without usage of any O-ring sealing which will release gas to affect vacuum purity when under high vacuum condition

**Vacuum Cavity:** V-Sorb unique *monolithic manifolds and pneumatic valve control system*, greatly reduce the dead volume; improve the adsorbate micro-change sensitivity; enhance pore size distribution analysis resolution; decrease connecting points; strengthen sealing performance and prolonged instrument life

**Measuring Ranges:** 0.01m<sup>2</sup>/g (Nitrogen)/ 0.005 m<sup>2</sup>/g (Krypton) to no known upper limit (surface area); 0.35 to 500nm (pore size)

**Accuracy:** repeatability errors  $\leq 1.0\%$

**Coolant Level Controller:** V-Sorb original *coolant level control system with temperature probe*, ensure the coolant level unchanged when compares with sample cells in the whole analysis process, completely eliminate the analysis errors caused by dead volume change

**Data Acquisition:** high-precision and high integration data acquisition and archiving system, minimal error, strong anti-interference ability

**Sample Ports:** two samples' analyzing and two samples' degassing concurrently

**Control System:** programmable pneumatic valve system with high integration and strong anti-interference ability, enhance instrument's stability and life

**Pressure Measurement:** imported *sectional measuring dual pressure transducer*, notably improve the measuring accuracy at low P/P<sub>0</sub> point, 0-1000Torr, 0-1Torr and 0-0.1Torr equipped.

**Transducer Accuracy:** imported *silicon thin film pressure transducer*, accuracy can reach 0.1% of real reading, better than 0.1% of F.S.(full-scale), far accurate than Pirani resistance vacuum gauge (general error is 10%-15%)

**Vacuum Pump:** built-in *bipolar Atlas Copco vacuum pump* controlled by patented software which can auto control pump's start/stop

**Dewar Vessel:** Japan brand *4 liter stainless steel Dewar* can provide at least 72 hours of unattended measurement without refilling the Dewar


**Partial Pressure:** P/P<sub>0</sub> controllable accuracy range is  $5 \times 10^{-8}$ -0.998

**Ultimate Vacuum:** bipolar pump:  $4 \times 10^{-2}$ Pa ( $3 \times 10^{-4}$ Torr), molecular pump:  $5 \times 10^{-8}$  mbar

**Sample Types:** powders, particle, fiber, flakes and other materials

**Adsorbate Gas:** high purity nitrogen ( $\geq 99.999\%$ ), Ar, Kr, CO, CO<sub>2</sub>, etc. non-corrosive gases are optional

**Data Reduction:** Windows®-based independent developed Gold APP Instruments software™, perfect versatility, produced full featured and multi-model reports



## Turbo Molecular Vacuum Pumped

### Surface Area & Porosity Analyzer **V-SORB 2800TP** Features

1. *Stainless steel micro welding vacuum manifolds* system can reduce the leak rate apparently, improve ultimate vacuum, without daily maintenance and enhance stability.
2. Micro welding vacuum system can decrease manifolds' volume, increases measuring sensitivity and raise accuracy.
3. Imported with *original packaging turbo molecular vacuum pump*, ultimate vacuum can reach  $5 \times 10^{-8}$ -0.998 mbar, maintenance free, convenient to use, reduce latter utilization cost.
4. Micro welding vacuum system, *high vacuum pneumatic valve* and high ultimate vacuum molecular pump unite perfectly, bring each accessories performance into full play, precision degree can reach until  $5 \times 10^{-8}$ .
5. Imported high vacuum pneumatic valve's sealing performance is up to  $5 \times 10^{-12}$  (Pa \*  $m^3/s$ ), can completely remove analyzed gas temperature change caused by valve heat.
6. Fully *metal-to-metal VCR* connecting can eliminate O-ring self gas release problem when under high vacuum circumstance, efficiently improve ultimate vacuum and analysis reliability.
7. Imported high-precision silicon thin film capacitive transducer, the core part, can realize auto zero by software; the *0.1 Torr full-scale (F.S.) transducer* assists micropore to achieve a more accuracy data for low P/P<sub>0</sub> points.
8. *Monolithic design* can configure as customer requests, benefits future functions extension and instrument maintenance.

9. Programmable pneumatic valve system with high integration and strong anti-interference ability, speed and easy for upgrade and maintenance.
10. *Fully auto and modularity design*, equipped with dedicated software, can realize unattended operation at night.
11. Multi calculating methods for data reduction provides all-round sample analysis options; powerful data archiving and searching system helps a lot for data management.



Japan  
Fujikin  
brand  
pneumatic  
valve



Inficon  
brand  
pressure  
transducers



# Testing Reports

Reports formats PDF+Excel+TXT+Word  
 Data cover BET/Langmuir/BJH/SF/T-plot/HK/DR/DA...  
 Selective P/Po points over 1000  
 Detailed adsorption & desorption points are available



## Surface Area and Pore Analysis Report


V-Sorb 2800TP Surface Area and Porosimetry Analyzer

### Sample Information

Sample ID: NO1	Model: V-Sorb 2800TP Surface Area and Poros
Submitter: micropore	Company: Gold APP Instruments China
Operator: Qiaoyun.Qin	Report Time: 2014-07-28

### Analysis Information

Mass: 0.13170 (g)	Pretreatment: 200 degree 2hrs
Method: Pore size	
Room Temp.: 25 degree	
Analysis Time:	



## Surface Area and Pore Analysis Report


V-Sorb 2800TP Surface Area and Porosimetry Analyzer

### Sample Information

Sample ID: NO1	Model: V-Sorb 2800TP Surface Area and Poros
Submitter: micropore	Company: Gold APP Instruments China
Operator: Qiaoyun.Qin	Report Time: 2014-07-28

### Analysis Information

Mass: 0.13170 (g)	Pretreatment: 200 degree 2hrs
Method: Pore size	
Room Temp.: 25 degree	
Analysis Time:	



www.app-one.com.cn

No.:

Manufacturer: APP Instrument

V-Sorb 2800TP Surface Area and Porosimetry Analyzer

Support: (8610)88099139




www.app-one.com.cn

No.:

Manufacturer: APP Instrument

V-Sorb 2800TP Surface Area and Porosimetry Analyzer

Support: (8610)88099139



# Accessory List

Sr. No.	Supplier	Parts	Qty.
1	Manufacturer Supplied	V-Sorb 2800TP Analyzer (with vacuum pump)	1 set
2		Rubber O-rings for Sample Cells Sealing	10 pcs
3		Spherical Sample Cells	10 pcs
4		V-shape Sample Funnel	10 pcs
5		Certificated Reference Material (large value)	10 g
6		Certificated Reference Material (medium value)	10 g
7		Certificated Reference Material (small value)	10 g
8		Copper Gas Pipe	2 m
9		Analysis Dewar (one 4L + one 520ml)	2 pcs
10		Fuse	2 pcs
11		RJ 45 Cable	2 pcs
12		Protective Gloves	1 pair
13		User Manual (English)	1 copy
14		Software CD (English)	1 copy
15		Po Cell	4 pcs
16		Filling Rod	8 pcs
17		Sample Cell Cleaning Brush	1 pc
18		Funnel Cleaning Brush	1 pc
19		Sample Weighting Cup	1 pc
20	Customer Prepared	Computer (Window® 2000/07/08/XP/Vista etc)	1 set
21		0.0001 precision balance	1 set
22		Printer (not a must)	1 set
23		Gas Regulator (should fit 1/8" NPT gas pipe, max reading is larger than 0.6Mpa.)	2 sets
24		Liquid Nitrogen (purity upper 99.999%)	1 pot
25		Power Cable	1 pc
26		He Gas (purity upper 99.999%)	1 cylinder
27		N <sub>2</sub> Gas (purity upper 99.999%)	1 cylinder

### **Headquarters**

Gold APP Instruments Corp. China  
Room 811, New Material Build,  
No.7<sup>th</sup>, Fenghui Mid. Rd.,  
Haidian District, Beijing 100094,  
P.R.China

Tel: 0086-10-88099138  
Mobile: 0086-18210009838  
Fax: +86-10-82132122  
Email: sales@jinaipu.com  
appone2008@hotmail.com

Skype:  
Gold-APP-Instruments  
We Chat & WhatsApp:  
0086-18210009838

### **Laboratory**

Room 601, New Material Building,  
No. 7th, Fenghui Mid. Rd.,  
Haidian Dist., Beijing 100094,  
P.R.China  
Tel: 0086-10-58711838  
Fax: 0086-10-58711838

### **Branch Offices**

Gold APP Instruments (Nanjing) Corp. China  
Room 512nd, No 4th Building,  
Mingfa Commerce Square, No. 99th,  
Yulan Rd., Yuhua District,  
Nanjing 210012,  
P.R.China  
Tel: 0086-25-58491095  
Fax: 0086-25-58491095

Gold iCON Instruments (Wuhan) Corp. China  
Room 5068, No. 1st Building,  
Huiyuan Block, No. 1st Rd.,  
Wuhan University Science Park,  
East Lake High-Tech Zone,  
Wuhan 430223,  
P.R.China  
Tel: 0086-27-59712850/1/2  
Fax: 0086-27-59712851 Ext.616

### **Webs**

[www.app-one.com.cn](http://www.app-one.com.cn)  
[www.jinaipu.com](http://www.jinaipu.com)



**Static volumetric principle analyzers:**  
[V-Sorb X800 series](#)

**Gas pycnometer true density analyzers:**  
[G-DenPyc X900 series](#)

**High pressure volumetric analyzers:**  
[H-Sorb X600 series](#)

**Sample preparation degassers:**  
[F-Sorb & V-Sorb](#)

*Our policy of continuous development may cause the information and specifications contained herein to change without notice or liability.*

### **A WORLD OF APPLICATION**

activated carbon, silica gel, active alumina oxide, molecular sieve, sepiolite, zeolite, alumina oxide, silicates, quartz, silicon carbide, lithium cobalt oxide, lithium manganese oxide, black lead, lithium nickel and cobalt, cobaltous oxide, lithium iron phosphate, lithium titanate, polymer, corrosion resister, silica, nano-calcium carbonate, zinc oxide, magnesium oxide, barium oxide, iron oxide, copper oxide, ferroferric oxide, ferrite, silver/iron/copper/ tungsten/nickel/aluminate powder, filler, inorganic filler, calcium carbonate, silica, deposited matter, suspended matter, titanium dioxide, rare earth, coal, cement, energy storage materials, catalyst, diatomaceous earth, cleansing agent, filter aid, superfine fiber, porous fabric, composite material, methane, coalbed gas etc.